

ELECTROSURGICAL GENERATORS

Force FX™ Electrosurgical Generator

with Instant Response™ technology



INSTANT RESPONSE™ TECHNOLOGY

Instant Response™ technology provides surgeons with improved performance at lower power settings, minimizing the risk of tissue damage and neuromuscular stimulation, and decreasing the need to “turn up the generator.”

A SMOOTHER CUT THROUGH ALL TISSUE TYPES

The Force FX™ generator adjusts automatically, responding to tissue changes, maintaining power delivery, and minimizing drag.

ADDED SAFETY AND RELIABILITY MINIMIZES ELECTROSURGICAL RISKS

Capacitive coupling is reduced by 30-50% when using Instant Response™ technology. This reduction is achieved by limiting the RMS voltage and the high frequency harmonics. Lower voltage means less neuromuscular stimulation and more precise delivery of energy to reduce collateral damage.



call 800.255.8522 to order

LIGASURE™ VESSEL SEALING SYSTEM
ELECTROSURGICAL GENERATORS
ELECTROSURGICAL ACCESSORIES
LAPAROSCOPIC INSTRUMENTS
SMOKE EVACUATION PRODUCTS
ARGON-ENHANCED ELECTROSURGERY
COOL-TIP™ RF SYSTEM
ULTRASONIC SURGICAL PRODUCTS
NEUROSURGICAL PRODUCTS
EDUCATION & TRAINING



What is Instant Response™ Technology?


Instant Response™ technology features an advanced feedback system that recognizes changes in tissue 200 times per second, and *adjusts* voltage and current accordingly to maintain appropriate power. This unique capability differentiates Instant Response™ generators from any other, as shown by their high power efficiency rating (PER). That's smart™ generator technology.

PER indicates the ability of a generator to accurately deliver the selected power into a wide range of tissue types. With a PER of 98, the Force FX™ generator provides a consistent surgical effect and has a higher PER than any other electrosurgical generator.*

Instant Response™ technology provides consistent power through different tissue types.

The Force FX™ generator provides surgeons with the advantage of accurately delivering the selected power setting, even through resistant tissue. The power curve below illustrates the Force FX™ generator's consistent output.

When a conventional generator encounters tissue changes, both power delivery and cutting efficiency are reduced. This means that higher power settings are required in order to cut through all types of tissue. With higher power settings, more energy is delivered, and tissue damage increases.



* PER values are approximate and are measured at 40 watts.

Unique Features of the Force FX™ Electrosurgical Generator

- Instant Response™ technology ensures that the power delivered remains virtually constant, regardless of the tissue type
- Improved performance at lower power settings minimizes the risk of tissue damage and neuromuscular stimulation
- Three internal microcontrollers reduce system reaction time and increase the system's processing speed
- Spray coagulation voltage of no more than 9000 volts peak-to-peak output for broad, but superficial coagulation with limited capacitive coupling
- A Power Efficiency Rating (PER) of approximately 98 for accurate and consistent cut performance
- Three cut modes, all controlled by Instant Response™ technology, offer surgeons a variety of choices
 - Low Cut for delicate tissue or laparoscopic cases
 - Pure Cut for a clean, precise cut
 - Blend for cutting with hemostasis
- Four coag modes
 - Desiccate for low voltage contact coagulation suitable in laparoscopic and delicate tissue work
 - Fulgurate (high crest factor) for efficient noncontact coagulation in most applications
 - Fulgurate (low crest factor) for lower voltage coagulation requirements
 - Spray for coagulating large tissue areas with superficial depth of necrosis
- Three bipolar modes
 - Precise, Standard, and Macrobipolar are controlled by the Instant Response™ system
 - Precise and Standard settings utilize low voltage to prevent sparking
- Versatile system that is uniquely compatible with other devices, including:
 - Force Argon™ II and Force GSU™ argon coagulation system
 - CUSA EXcel™ and CUSA™ 200 ultrasonic surgical aspirators
 - OptiMumm™ smoke evacuator, through a direct cable link
 - Valleylab VLCM bipolar current monitor
- Compatible with and used as the electrosurgical energy source for:
 - Dyonics* Control RF arthroscopic ablation system
 - Dyonics* Electroblade™ rotary resection system
 - Cook Vascular Perfecta™ EDS pacemaker lead extraction system
- Compatible with, and the exclusive electrosurgical generator for, the Computer Motion Hermes™ Voice Command System

Force FX™ Electrosurgical Generator

Technical Specifications (110-120V)

OUTPUT WAVEFORMS

Bipolar
 Precise: 470 kHz sinusoid
 Standard: 470 kHz sinusoid
 Macro: 470 kHz sinusoid

Monopolar Cut

Low: 390 kHz sinusoid. Similar to the Pure Cut mode except the maximum voltage is limited to a lower value.
Pure: 390 kHz sinusoid
Blend: 390 kHz bursts of sinusoid, recurring at 27 kHz intervals. 50% duty cycle envelope.

Monopolar Coag

Desiccate: 240 kHz sinusoid repeated at 39 kHz. 8% duty cycle.
Fulgurate: 390 kHz damped sinusoidal bursts with a repetition frequency of 30 or 57 kHz into 500 ohms
Spray: 390 kHz damped sinusoidal bursts with a randomized repetition centered at 28 kHz. Frequencies include 21 kHz < f < 35 kHz. Output is further modulated by a random 250 Hz envelope with a variable duty cycle.

Output power changes by less than 15% or 5 watts, whichever is greater, as the line voltage varies from 104-132 volts and 208-264 volts (at rated load).

LOW FREQUENCY LEAKAGE (50-60 HZ)

Source current, patient leads, all outputs tied together.

- Normal polarity, intact chassis ground < 10 µA
- Normal polarity, ground open < 50 µA
- Reverse polarity, ground open < 50 µA
- Sink current, 140V applied, all inputs < 50 µA

PER = 98

PER (Power Efficiency Rating) is the measure of an electrosurgical generator's ability to accurately deliver the selected power into a wide range of tissue types.

WEIGHT AND DIMENSIONS

Height: 11.1 cm (4-3/8 in.)
Width: 35.6 cm (14 in.)
Length: 43.9 cm (17 in.)
Weight: < 8.1 kg (< 18 lbs)

INPUT POWER REQUIREMENTS

Operating range is 85 to 132 AC volts. Maximum current is 7 amperes in Cut and 4 amperes in Coag.

HIGH FREQUENCY LEAKAGE

Bipolar: Less than 60 mA_{rms}
Monopolar: Less than 150 mA_{rms}

REM™ CONTACT QUALITY MONITORING SYSTEM

Measurement Frequency: 80 kHz ± 10 kHz
Measurement Current: Less than 10 µA

Acceptable Resistance Ranges:
 REM™ pad – 5-135 ohms
 Non-REM™ pad – less than 20 ohms

Acceptance range is 5-135 ohms after REM PolyHesive™ II return electrode is applied.

Adaptive REM™

REM™ trip is baseline impedance plus 40%. For example, if the baseline impedance is 30 ohms, the upper level trip is approximately 42 ohms. If the pad-patient impedance falls below the baseline impedance, a new baseline is established.

Meets UL and cUL specifications.



OUTPUT CHARACTERISTICS

	MODE	MAXIMUM P-P VOLTAGE	RATED LOAD (OHMS)	MAXIMUM POWER (WATTS)	CREST FACTOR* (TYPICAL)
Bipolar	Precise	450	100	70	1.5
	Standard	320	100	70	1.5
	Macro	750	100	70	1.5
Monopolar Cut	Low	1350	300	300	1.5
	Pure	2300	300	300	1.5
	Blend	3300	300	200	2.5
Monopolar Coag	Desiccate	3500	500	120	5
	Fulg. High Crest Factor	8500	500	120	7.0
	Fulg. Low Crest Factor	6900	500	120	5.5
	Spray	9000	500	120	8
CEM™ Mode	Monopolar Cut (Low)	1000	300	100	1.5
	Monopolar Coag (Desiccate)	3500	500	70	5

*Crest Factor is an indicator of a waveform's ability to coagulate without cutting.

ORDER INFORMATION

CATALOG NUMBER	DESCRIPTION	ORDER QUANTITY
Force FX-C	Microcontroller-based isolated electrosurgical generator, designed for all general surgical procedures. Unit includes the Valleylab adaptive REM™ system and Instant Response™ technology.	1 each

Valleylab
 a division of Tyco Healthcare Group LP
 5920 Longbow Drive
 Boulder, CO 80301-3299 USA
 800 255 8522

